



**NSAI**  
Standards

Standard Recommendation  
S.R. CEN/TS 15634-2:2012

# Foodstuffs - Detection of food allergens by molecular biological methods - Part 2: Celery (*Apium graveolens*) - Qualitative determination of a specific DNA sequence in cooked sausages by real-time PCR

## S.R. CEN/TS 15634-2:2012

*Incorporating amendments/corrigenda/National Annexes issued since publication:*

The National Standards Authority of Ireland (NSAI) produces the following categories of formal documents:

I.S. xxx: Irish Standard – national specification based on the consensus of an expert panel and subject to public consultation.

S.R. xxx: Standard Recommendation - recommendation based on the consensus of an expert panel and subject to public consultation.

SWIFT xxx: A rapidly developed recommendatory document based on the consensus of the participants of an NSAI workshop.

*This document replaces:*

*This document is based on:*  
CEN/TS 15634-2:2012

*Published:*  
10 February, 2012

This document was published  
under the authority of the NSAI  
and comes into effect on:  
10 February, 2012

**ICS number:**

07.100.30  
67.050

**NSAI**  
1 Swift Square,  
Northwood, Santry  
Dublin 9

T +353 1 807 3800  
F +353 1 807 3838  
E standards@nsai.ie  
W NSAI.ie

**Sales:**  
T +353 1 857 6730  
F +353 1 857 6729  
W standards.ie

Údarás um Chaighdeáin Náisiúnta na hÉireann

ICS 07.100.30; 67.050

English Version

**Foodstuffs - Detection of food allergens by molecular biological methods - Part 2: Celery (*Apium graveolens*) - Qualitative determination of a specific DNA sequence in cooked sausages by real-time PCR**

Produits alimentaires - Détection des allergènes alimentaires par des méthodes d'analyse de biologie moléculaire - Partie 2: Céleri (*Apium graveolens*) - Détermination qualitative d'une séquence d'ADN spécifique dans des saucisses cuites par PCR en temps réel

Lebensmittel - Nachweis von Lebensmittelallergenen mit molekularbiologischen Verfahren - Teil 2: Sellerie (*Apium graveolens*) - Qualitative Bestimmung einer spezifischen DNA-Sequenz in Brühwürsten mittels Real-time-PCR

This Technical Specification (CEN/TS) was approved by CEN on 15 November 2011 for provisional application.

The period of validity of this CEN/TS is limited initially to three years. After two years the members of CEN will be requested to submit their comments, particularly on the question whether the CEN/TS can be converted into a European Standard.

CEN members are required to announce the existence of this CEN/TS in the same way as for an EN and to make the CEN/TS available promptly at national level in an appropriate form. It is permissible to keep conflicting national standards in force (in parallel to the CEN/TS) until the final decision about the possible conversion of the CEN/TS into an EN is reached.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.



EUROPEAN COMMITTEE FOR STANDARDIZATION  
COMITÉ EUROPÉEN DE NORMALISATION  
EUROPÄISCHES KOMITEE FÜR NORMUNG

**Management Centre: Avenue Marnix 17, B-1000 Brussels**

<b>Contents</b>		<b>Page</b>
<b>Foreword.....</b>		<b>3</b>
<b>1</b>	<b>Scope .....</b>	<b>4</b>
<b>2</b>	<b>Principle .....</b>	<b>4</b>
<b>3</b>	<b>Reagents .....</b>	<b>4</b>
<b>4</b>	<b>Apparatus and equipment .....</b>	<b>6</b>
<b>5</b>	<b>Analysis steps .....</b>	<b>6</b>
<b>6</b>	<b>Validation status and performance criteria .....</b>	<b>11</b>
<b>7</b>	<b>Sample type and amounts .....</b>	<b>14</b>
<b>8</b>	<b>Limit of detection .....</b>	<b>14</b>
<b>9</b>	<b>Interferences .....</b>	<b>14</b>
<b>10</b>	<b>Test report .....</b>	<b>14</b>
<b>11</b>	<b>Method performance studies .....</b>	<b>15</b>
<b>Bibliography .....</b>		<b>17</b>

## **Foreword**

This document (CEN/TS 15634-2:2012) has been prepared by Technical Committee CEN/TC 275 “Food analysis - Horizontal methods”, the secretariat of which is held by DIN.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CEN [and/or CENELEC] shall not be held responsible for identifying any or all such patent rights.

According to the CEN/CENELEC Internal Regulations, the national standards organizations of the following countries are bound to announce this Technical Specification: Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom.

## 1 Scope

This Technical Specification specifies a method for the qualitative detection of celery (*Apium graveolens*) in emulsion-type sausages (e.g. Frankfurter, Wiener).

Real-time PCR detection of celery is based on an 101 bp (base pair) sequence from the gene of the mannitol dehydrogenase (GenBank Acc. No. AF067082) of celery (*Apium graveolens*).

The method has been validated on emulsion-type sausages (Bavarian "Leberkäse") spiked with celery. For this purpose meat batter containing mass fractions of 50 % pork meat, 25 % pork fat, 23 % crushed ice and 1,8 % of a mixture of sodium chloride, nitrite, nitrate, phosphates and ascorbates was prepared according to a standard procedure for emulsion-type sausage. The meat batter was spiked with either ground celery seeds or celery root powder to 1000 mg/kg. Lower spiking levels were obtained by diluting with celery-free meat batter. The batter was stuffed into casings and heated at 65 °C for 60 min [2].

## 2 Principle

Total DNA from emulsion-type sausages are isolated from the sample matrix. DNA is released from the sample matrix using the cetyltrimethylammonium bromide (CTAB) approach. Potential PCR inhibitors are removed from the isolated DNA by purification with solid phase columns. Real-time PCR is used to detect, amplify and quantify a celery specific sequence. The real time PCR method involves a fluorescence approach with a sequence specific hydrolysis probe [1], [2].

## 3 Reagents

### 3.1 General

The following general conditions for analysis shall be followed, unless specified differently. Use only analytical grade reagents suitable for molecular biology. Reagents shall be stored in small aliquots to minimise the risk of contamination. All water shall be free from DNA and nucleases, e.g., double distilled or equivalent (molecular grade). Solutions shall be prepared by dissolving the appropriate reagents in water and autoclaving, unless specified differently.

### 3.2 Extraction reagents

**3.2.1 Chloroform**, CAS 66-67-3.

**3.2.2 Ethanol**, volume fraction  $\varphi$  = 70 %, CAS 64-17-5.

**3.2.3 Ethylenediaminetetraacetic acid disodium salt (Na<sub>2</sub>EDTA)**, CAS 6381-92-6.

**3.2.4 Cetyltrimethylammoniumbromide (CTAB)**, CAS 57-09-0.

**3.2.5 Hydrochloric acid**,  $\varphi$  = 37 %, CAS 7647-01-0.

**3.2.6 Isoamyl alcohol**, CAS 123-51-3.

**3.2.7 Isopropanol**, CAS 67-63-0.

**3.2.8 Proteinase K**, EC 3.4.21.64.

**3.2.9 Sodium chloride**, CAS 7647-14-5.

This is a free preview. Purchase the entire publication at the link below:

[Product Page](#)

- 
- Looking for additional Standards? Visit Intertek Inform Infostore
  - Learn about LexConnect, All Jurisdictions, Standards referenced in Australian legislation
-